



Presence of blood or semen in a qfFN sample can sometimes lead to a falsely elevated result, however results under the specified threshold for treatment at your facility can still be considered valid in either of these situations.

Before undertaking a qfFN swab – think – “will the result of this change my management plan?”

*Manipulation of the cervix (from a vaginal examination, swab or vaginal ultrasound probe etc.) may interfere with your qfFN sample. It therefore must be the first clinical procedure you do to ensure you have a reliable result.



Use the QUIPP App

Search QUIPP for Android or App store, or go to www.quipp.org